



[www.smaio.com/ifu](http://www.smaio.com/ifu)



This document is valid only on the date printed. If unsure of the print date, please re-print to ensure use of the latest revision of the IFU available on the website [www.smaio.com/ifu](http://www.smaio.com/ifu).

Instructions for use only for the United States of America (USA).

Caution: Federal Law (USA) restricts this device to sale by or on the order of a licensed health practitioner

### 1. INTENDED USE

This instrument, herein referred to as “device”, is part of a kit dedicated to placement of a specific product in the KHEIRON range of implants. This device is re-usable, re-sterilizable and must be handled by healthcare professionals trained in its conditions for use. It must be handled with care at all stages of use, storage and maintenance.

### 2. INDICATIONS

This device belongs to the KHEIRON product range defined by a nomenclature described in the surgical technique\* provided with the implants and/or instruments of the product range. These devices can be used only with the KHEIRON range of implants.

### 3. CONTRAINDICATIONS

- Do not use this device if it is no longer functional.
- Do not use this device outside of its planned scope of application.
- Allergies, intolerances and/or hypersensitivity to the component materials of the instrument (see the KHEIRON surgical technique).

### 4. STORAGE CONDITIONS

Store the devices in a clean place, at room temperature in their original packaging or in the box provided for this purpose by SMAIO.

### 5. PRECAUTIONS AND WARNINGS

- Delivered instruments are NON-STERILE and are intended to be cleaned and sterilized according to the instructions given in §6. **Knowledge and experience in spinal surgery are prerequisites.**
- These instruments do not present any known risk of interference with other medical equipment.

### 6. CLEANING / STERILIZATION

- Whether they come directly from their original packaging or from the tray of use, each device must be cleaned and decontaminated in conformity with legislation in force prior to sterilization.
- Before use, these products must be sterilized by steam autoclaving in appropriate containers in accordance with ISO standard ISO 17665-1.
- SMAIO recommends cleaning the non-sterile devices by combined manual and automatic cleaning, using a heat disinfectant which complies with standard EN ISO 15883-1, used with an alkaline cleaning product (pH ≤ 10) as per the validated method described in the tables below.
- The aforementioned parameters are valid only for devices in the product range sterilized in the corresponding box. Any other configuration used would jeopardize the manufacturer's validation.
- The instructions provided below have been validated by SMAIO as being capable of preparing a medical device for reuse. The processor remains responsible for ensuring that the processing, as actually performed using equipment, materials and personnel at the processing facility, achieves the desired result. This requires verification and/or validation and routine process monitoring.
- Clean the device until there is no visual contamination of the instruments directly after application (within a maximum of 2 h of use).

<b>Manufacturer:</b> SMAIO <b>Method:</b> Steam autoclaving <b>Devices:</b> KHEIRON Spinal Fixation System	
<b>WARNINGS</b>	<ul style="list-style-type: none"> <li>Devices made of stainless steel, silicone, and POM (Polyacetal) must NOT BE processed with chlorinated products but can safely be immersed in NaOH (NaOH – 1 mol/L) for 1h.</li> <li>Devices made of radel or propylux can safely be immersed in NaOH or sodium hypochlorite.</li> <li>In all cases, do not use a wire brush or abrasive, and handle these devices with gloves throughout the various processes and uses, during which they should be placed on trays suitable for the various cleaning and decontamination steps.</li> </ul>
<b>Limitations on Processing</b>	Non applicable

INSTRUCTIONS																													
<b>Initial treatment at point use</b>	Place the instruments on a tray or in a suitable basket for the cleaning and decontamination steps.																												
<b>Preparation before cleaning</b>	<p>Disassemble the devices when possible (see disassembly instructions). Open the devices, when possible.</p> <p>Fully immerse the devices in an enzymatic solution for the time recommended by the cleaning agent manufacturer.</p> <p>Validated parameters are: “Neodisher® MediClean Forte” 0.5% cleaning agent; room temperature [20 °C; 25 °C] [68 °F; 77 °F] ; for 10 minutes.</p>																												
<b>Pre-cleaning: manual</b>	<p>Remove any visible residues with a soft brush (non-wire) or a long neck brush capable of reaching the cavities. Brush for at least 2 minutes. Carefully clean each cavity.</p> <p>Remove the device from the enzymatic solution. Abundantly rinse the devices with tap water, at room temperature [20°C; 25°C] [68 °F; 77 °F] for 1 minute and repeat the operation at least 3 times. Carefully rinse each lumen, rough or hard-to-reach area, cavity, etc.</p>																												
<b>Cleaning: Automated</b>	<p>Place the instruments in an opened position in a suitable basket for the cleaning and decontamination step.</p> <p><b>Parameters validated with the cleaning/disinfectant HAMO ECOLINE LM-25 and the “VARIO TD LIQUID” program</b></p> <table border="1" style="width: 100%; text-align: center;"> <thead> <tr> <th>Steps</th> <th>Time (min)</th> <th>Temperature</th> <th>Cleaning agent</th> </tr> </thead> <tbody> <tr> <td>Pre-cleaning</td> <td>2 :00</td> <td>&lt; 45 °C &lt; 113 °F</td> <td>Tap water</td> </tr> <tr> <td>Cleaning</td> <td>5 :00</td> <td>55 °C - 131°F</td> <td>Neodisher® MediClean Forte 2mL per litre</td> </tr> <tr> <td>Neutralisation</td> <td>2:00</td> <td>&lt; 45 °C &lt; 113 °F</td> <td>Tap water</td> </tr> <tr> <td>Rinsing</td> <td>2:00</td> <td>&lt; 45 °C &lt; 113 °F</td> <td>Tap water</td> </tr> <tr> <td>Disinfection</td> <td>5:00</td> <td>90 °C - 194°F</td> <td>Osmosis water</td> </tr> <tr> <td>Drying</td> <td>22:00</td> <td>80 °C - 176°F</td> <td>N/A</td> </tr> </tbody> </table>	Steps	Time (min)	Temperature	Cleaning agent	Pre-cleaning	2 :00	< 45 °C < 113 °F	Tap water	Cleaning	5 :00	55 °C - 131°F	Neodisher® MediClean Forte 2mL per litre	Neutralisation	2:00	< 45 °C < 113 °F	Tap water	Rinsing	2:00	< 45 °C < 113 °F	Tap water	Disinfection	5:00	90 °C - 194°F	Osmosis water	Drying	22:00	80 °C - 176°F	N/A
Steps	Time (min)	Temperature	Cleaning agent																										
Pre-cleaning	2 :00	< 45 °C < 113 °F	Tap water																										
Cleaning	5 :00	55 °C - 131°F	Neodisher® MediClean Forte 2mL per litre																										
Neutralisation	2:00	< 45 °C < 113 °F	Tap water																										
Rinsing	2:00	< 45 °C < 113 °F	Tap water																										
Disinfection	5:00	90 °C - 194°F	Osmosis water																										
Drying	22:00	80 °C - 176°F	N/A																										
<b>Maintenance, inspection and Testing</b>	After the cleaning cycle, check each device for any visible oil residue. If residue is found, repeat each cleaning step until there is no visible residue left.																												
<b>Packaging</b>	The devices are placed in the SMAIO kit. The kit is packaged in a double Pasteur wrap.																												

<b>Sterilization</b>	Assembled devices must NOT be re-assembled before the sterilization cycle. The sterilization cycles are the following, depending on the geographic area:			
		<b>European Union</b>	<b>United Kingdom</b>	<b>United States</b>
	<b>Method: steam</b>	Yes	Yes	Yes
	<b>Pre-vacuum</b>	Yes	Yes	Yes
	<b>Temperature</b>	134 °C 273.2 °F	134 °C 273.2 °F	132 °C 269.6 °F
	<b>Time</b>	18 min	3 min	4 min
	<b>Drying time</b>	45 min	45 min	45 min
	Wrap tray with a towel placed between tray and FDA cleared wrap.			
<b>Storage</b>	Not applicable			
<b>Additional information</b>	Look for any signs of premature instrument wear after sterilization. If such is the case, DO NOT use them and inform SMAIO (see §11).			
<b>Manufacturer contact</b>	See §11			

## 7. INSTRUCTIONS FOR USE

Instructions for placement of implants are detailed in the surgical technique.\*

If in doubt concerning use, cleaning, decontamination, sterilization or discarding of an instrument, do not hesitate to contact the SMAIO customer service department or its qualified distributor.

### a) Before the operation

- Carefully read the instructions for use of the KHEIRON implants to be used in association with the KHEIRON instruments, along with the surgical technique\*.
- Carefully read the documents provided with the instrumentation kit.
- Prepare all implants and instruments required for the procedure and ensure that nothing is missing.
- Check their cleanliness, integrity and functionality.
- Only SMAIO implant placement instruments, designed and provided by SMAIO, should be used with the implants.
- Before a first implantation, the surgeon and their surgical team should practice handling the instruments to become familiar with the equipment.

### b) During the operation

- The procedure should be performed by a surgeon who has received the necessary spinal surgery training.
- Observe the different phases described in the surgical technique\*.
- Do not use this device for purposes other than those for which it is intended (see the surgical technique\*)

- The implantation must be performed solely with the instruments provided for this purpose and according to the indications of the surgical technique\*.

### c) After the operation

- All materials must be cleaned and decontaminated immediately before transport (see §6)
- Discard any products which have not functioned properly and ask for a replacement
- Disassemble this product if necessary according to the assembly/disassembly instructions\*, scrupulously clean and decontaminate it, then sterilize it before returning it or putting it back into service (see §6).

## 8. PRODUCT DISPOSAL

To discard a product following an error in product storage or use, the instrument must follow the pathway for removal of hospital waste products in compliance with procedures in force within the institution.

## 9. GUARANTEE

In the event of a defect, contact the SMAIO customer service department and return the defective device, cleaned, decontaminated and together with the liaison form available on SMAIO website.

## 10. VIGILANCE

Any adverse events linked to the use of the device should be reported to SMAIO at [usvigilance@smaio.com](mailto:usvigilance@smaio.com) and to the FDA.

## 11. MANUFACTURER CONTACT



S.M.A.I.O  
2, Place Berthe Morisot – Parc Technologique  
69800 SAINT-PRIEST – France  
Tel: +33 (0)4 69 84 23 02  
Website: [www.smaio.com](http://www.smaio.com)

\* Documentation available on request on the [www.smaio.com](http://www.smaio.com) website.

## 12. GLOSSARY OF SYMBOLS USED IN ALL KHEIRON LABELING

Symbol	Meaning of symbol	Reference
	Medical Device	ISO 15223-1, 5.7.7
	Do not re-use	ISO 15223-1, 5.4.2
	Non-sterile	ISO 15223-1, 5.2.7
	Caution	ISO 15223-1, 5.4.4
	Do not use if package is damaged	ISO 15223-1, 5.2.8
	Keep away from sunlight	ISO 15223-1, 5.3.2
	Keep dry	ISO 15223-1, 5.3.4
	Consult instructions for use	ISO 15223-1, 5.4.3
	Manufacturer	ISO 15223-1, 5.1.1
	CE Marking of conformity	MDD 93/42/CE Annex XII
	Date of manufacture	ISO 15223-1, 5.1.43
	Catalogue number	ISO 15223-1, 5.1.46
	Batch code	ISO 15223-1, 5.1.54
	Serial number	ISO 15223-1, 5.1.7
	Federal Law (USA) restricts this device to sale by or on the order of a physician	21 CFR 801.109